

SORBALGON ALGINATE DRESSING



MAIN SPECIFICATIONS

- Sorbalgon polysorbate is made of calcium alginate fibers and forming a non-adherent hydrophilic gel provided with the sodium salt which is in the blood and wound exudate. The gel forming property of promoting wound healing process.
- The gel with hemostatic properties for bleeding wounds fills the wound cavity and creates favorable conditions for healing, removing bacteria from the wound and keeping them within the gel. Due to the physical and chemical properties pad has a high absorption capacity which is very important for the treatment of moderate to high exuding wounds. Adapts very well to contours and even some forms of difficult wounds. Under surveillance can remain in the wound a few days and should be changed when it is completely transformed into gel.
- Primarily for the treatment of external wounds of all kinds, but particularly suitable for the treatment of wounds with blood and exudate, gel forming properties favor the process of wound healing, for example in the venous ulcers, bedsores, abscesses, boils, burns, wounds difficult to heal caused by accident or tumors surgery.
- Especially suitable for treatment of bleeding or exuding wounds, since the gel formation properties conducive to wound healing process.
- In the case of deep wounds the dressing dens or calcium alginate may be used as the buffer in places difficult to access. Through absorption of wound exudate the dressing quickly becomes a gel which ensures contact with the wound bed.

- The dressings in the form of sterile swab, are available in three different sizes to basic needs and using economic and available in the form of buffers in strips 30 cm long, especially for narrow fistula passage points and channels.

Code	Description	Selling unit
01046037	SORBALGON ALGINATE DRESSING 5X5CM	CX. 10
01046017	SORGALGON ALGINATE DRESSING 10X10 CM	CX. 10
01046042	SORGALGON ALGINATE DRESSING 10X20 CM	CX. 5
9995901 (H)	SORBALGON T ALGINATE DRESSING 1G/30 CM	CX. 3
01046065	SORBALGON T ALGINATE DRESSING 2G/30 CM	CX. 5