

HYDROSORB GEL HYDROGEL



MAIN SPECIFICATIONS

- For residue-free moist wound treatment.
- 60% water content - immediately provides a moist environment to promote wound healing.
- High biocompatibility, conformable fit.
- Permanently moist so does not adhere to the wound.
- Improves patient comfort due to the natural anaesthetic and cooling effect of water on the skin.
- Imprinted 1cm square grid.
- Aids exact wound documentation.
- Semi-permeable, polyurethane film cover.
- Protects wound from penetration of moisture and bacteria.
- Allows free passage of gases and moisture vapour.
- Leaves wound residue-free for immediate assessment.
- Sustained absorbency and transparency.
- Exudate and malodours remain contained

Code	Description	Selling unit
01046045	HYDROSORB GEL 15 GR	CX. 10