

HYDROFILM PLUS DRESSING



MAIN SPECIFICATIONS

- Reliable transparent adhesive film dressings that protect wounds from bacteria and water.
- Breathable film.
- Hypoallergenic adhesive.
- Transparent, waterproof.
- Conform to body contours with a low profile.
- Numbered and colored application guides.
- Hydrofilm Plus includes an absorbent, non-stick pad.
- Prevents moisture build-up and dressing peel-off.
- Secures comfortably on skin for longer wear time.
- Allows wound assessment and showering without dressing removal.
- Secure and reliable fit.
- Easy-to-follow steps that can be performed with one gloved hand.
- Hydrofilm Plus wicks exudate away without sticking to wound.
- Hydrofilm should not be used as a primary dressing on clinically infected, bleeding or heavily secreting wounds.

- Hydrofilm is for use as a post-operative and primary dressing to protect against secondary infection and to cover non-exudating wounds.
- Hydrofilm is also ideal for use as a secondary dressing for fixation of catheters or cannulae, or as a secondary dressing to cover primary calcium alginates.
- Hydrofilm Plus is for use in the care and protection of light to moderate exudating wounds including minor cuts, burns and scrapes.

| Code | Description | Selling unit |
|-------------|---|--------------|
| 01046050 | HYDROFILM PLUS DRESSING WITH PAD 9X10 CM | CX. 50 |
| 01046034 | HYDROFILM PLUS DRESSING WITH PAD 9X15 CM | CX. 25 |
| 01046062 | HYDROFILM PLUS DRESSING WITH PAD 5X7,2 CM | CX. 50 |
| 01046063 | HYDROFILM PLUS DRESSING WITH PAD 10X12 CM | CX. 25 |
| 01046061 | HYDROFILM PLUS DRESSING WITH PAD 10X20 CM | CX. 25 |
| 6857790 (H) | HYDROFILM PLUS DRESSING WITH PAD 10X25 CM | CX. 25 |
| 6857800 (H) | HYDROFILM PLUS DRESSING WITH PAD 10X35 CM | CX. 25 |