

# HYDROFILM DRESSING



## MAIN SPECIFICATIONS

- Hydrofilm is a self adhesive post operative dressing which provides reliable protection of epithelialisation and fixation of medical devices.
- Hydrofilm Plus is a self adhesive post operative dressing, with a highly absorbent wound pad, to protect exuding wounds.
- Waterproof and impermeable to bacteria.
- High MVTR for breathability.
- Transparent.
- Suitable for sensitive skin.
- Easy to apply.

| Code     | Description                               | Selling unit |
|----------|---|--------------|
| 01046028 | HYDROFILM DRESSING WITHOUT PAD 6X7 CM     | CX. 100      |
| 01046030 | HYDROFILM DRESSING WITHOUT PAD 10X12,5 CM | CX. 100      |
| 01046046 | HYDROFILM DRESSING WITHOUT PAD 10X15 CM   | CX. 50       |