

# HYDROTUL HYDROACTIVE DRESSING



## MAIN SPECIFICATIONS

- Hydrotul is a hydroactive, non-occlusive tulle dressing, impregnated with a non-medicated triglyceride ointment, which contains carboxy-methol cellulose (CMC) for moist wound healing.
- When in contact with wound exudate, the hydrocolloid particles embedded in the ointment mass absorb and store wound moisture forming a gel at the wound surface. This gel provides an optimal moist environment for managing acute and chronic wounds.
- Hydrotul is indicated for wounds in the granulation and epithelialisation phases of wound healing, which require a moist wound environment.
- The dressing can be left on the wound for several days and may be covered with a secondary dressing such as Zetuvit or PermaFoam to absorb wound exudate.

Code	Description	Selling unit
01046048	HYDROTUL HYDROACTIVE DRESSING 5X5 CM	CX. 10
01046049	HYDROTUL HYDROACTIVE DRESSING 10X12 CM	CX. 10
4995841 (H)	HYDROTUL HYDROACTIVE DRESSING 15X20 CM	CX. 10
4995851 (H)	HYDROTUL HYDROACTIVE DRESSING 10X12 CM	CX. 50