

ATRAUMAN SILVER DRESSING



MAIN SPECIFICATIONS

- For atraumatic wound care.
- Hydrophobic polyester fibres with 1mm pore size.
- Prevents granulation tissue from penetrating the dressing, minimising pain and trauma on removal.
- Allows the flow of exudate and air.
- Provides high conformability.
- Neutral triglyceride impregnation.
- Does not contain Vaseline or Paraffin, leaving the wound residue free.
- Keeps wound margin soft and supple, avoiding scar tissue contraction.
- Enhances non-adherent properties, minimising pain and trauma on removal.
- “Touch & Kill” active antibacterial effect.
- When in contact with wound exudate, Atrauman Ag forms silver ions on its surface, killing bacteria on contact.
- Wound exudate, dead bacteria and endotoxins are absorbed by the secondary dressing.
- Silver ions are kept within the dressing where they bind to and destroy bacteria. Due to this local restriction of silver ion release, cytotoxicity of Atrauman Ag is very low.

Code	Description	Selling unit
01046018	ATRAUMAN SILVER DRESSING 5X5 CM	CX. 10
01046019	ATRAUMAN SILVER DRESSING 10X10 CM	CX. 10
4995369 (H)	ATRAUMAN SILVER DRESSING 10X20 CM	CX. 30