

# COSMOPOR DRESSING POSTOPERATIVE



## MAIN SPECIFICATIONS

- For sterile dressing of post operative wounds and/or minor injuries.
- Viscose Core that provides high absorbency and cushioning.
- Polyethylene net contact layer that prevents adherence of the dressing to the wound.
- Hypoallergenic fibres, permeable to air and water, that prevents irritation of surrounding skin.
- Enables the natural functions of the skin to continue.

Code	Description	Selling unit
01046033	COSMOPOR E DRESSING POSTOPERATIVE 6X10 CM	CX. 25
01046040	COSMOPOR E DRESSING POSTOPERATIVE 6X15 CM	CX. 25
01046041	COSMOPOR E DRESSING POSTOPERATIVE 10X35 CM	CX. 25
01046025	COSMOPOR E DRESSING POSTOPERATIVE 8X10 CM	CX. 25
01046026	COSMOPOR E DRESSING POSTOPERATIVE 8X15 CM	CX. 25

01046027	COSMOPOR E DRESSING POSTOPERATIVE 8X20 CM	CX. 25
01046014	COSMOPOR E DRESSING POSTOPERATIVE 7,2X5 CM	CX. 50
01046000	COSMOPOR E DRESSING POSTOPERATIVE 10X20 CM	CX. 25
9008775 (H)	COSMOPOR E DRESSING POSTOPERATIVE 10X25 CM	CX. 25