

COSMOPOR STRIP DRESSING



MAIN SPECIFICATIONS

- The wound dressing with the special open-cell nonwoven.
- Cosmopor Strip is a hypoallergenic adhesive plaster made of soft fleece fabric with synthetic rubber adhesive and a good absorbent dressing. Due to the elasticity of the nonwoven, Cosmopor Strip can be well adapted to conical body parts. It is ideal for the treatment of minor injuries in patients with very sensitive skin.
- For the treatment of minor injuries, particularly on joints and conical parts of the body; in patients with sensitive skin or on particularly sensitive body areas. For the application the skin should be dry and clean.

Code	Description	Selling unit
01046005	COSMOPORE STRIP DRESSING 5 CM X4 M	CX. 1
01046006	COSMOPORE STRIP DRESSING 10X4 CM	CX. 10
01046007	COSMOPOR STRIP DRESSING 10X8 CM	CX. 10