

# STERI-STRIP SUTURE / CLOSURE STRIP



## MAIN SPECIFICATIONS

- Steri-Strip is a sterile micropore adhesive with polyester filaments, which promotes skin closing quickly and easily.
- It is pressure sensitive, has good initial adhesion, providing great savings in time compared to conventional sutures.
- Being a non-invasive adhesive reduces the likelihood of trauma to the tissue and the consequent formation of dermal edema.
- It provides faster healing and reduces the risk of infections caused by invasive sutures - cutaneous sutures packaged individually.
- With high adhesion.
- It allows appropriate sweating skin.
- Adhesive skin sutures, sterile.
- Back viscose rayon nonwoven.
- Reinforced strips of polyester.
- Hypoallergenic acrylic adhesive.
- Latex-free.
- Indications: closing surgical incisions; Approach uninfected wounds edges; Closing clean superficial

sutures (abrasions, lacerations); As a complement to conventional sutures, increasing the tensile strength and support to the wound after removing the conventional suture.

- Available in various sizes.

Code	Description	Selling unit
01091011	STERI STRIP 3X75 MM	CX. 50
01091012	STERI STRIP 6X75 MM	CX. 50
01091013	STERI STRIP 6X38 MM	CX. 50
01091014	STERI STRIP 6X100 MM	CX. 50
01091015	STERI STRIP 12X100 MM	CX. 50
01091017	STERI STRIP 3X75 MM	CX. 5
01091018	STERI STRIP 6X75 MM	CX. 3