

# TEGADERM PAD DRESSING WITH NON-ADHERENT PAD



## MAIN SPECIFICATIONS

- Convenient, all-in-one sterile dressing reduces application time.
- Waterproof barrier which allows patients to bathe or shower without removing dressing.
- Unique, non-adherent pad won't stick to the wound, allowing normal wound healing processes to occur with less pain and trauma.
- Breathable film allows skin to function normally with good exchange of moisture vapor and oxygen.
- Conformable dressing is easy to apply, and stretches with the skin to promote longer wear time.
- This island dressing is a sterile, waterproof, viral and bacterial barrier which consists of a non-adherent absorbent pad bonded to a larger thin film transparent dressing (In vitro testing shows that the transparent film provides a viral barrier from virus 27 nm in diameter or larger while the dressing remains intact without leakage).
- Post-operative wounds; Acute wounds such as abrasions and lacerations; Light-to-moderately-draining chronic wounds; Superficial and partial-thickness burns; I.V. catheter sites (peripheral and central lines); To secure or cover devices; Secondary dressing applications: Over hydrogels, alginates, and gauze

Code	Description	Selling unit
------	-------------	--------------

01091016	TEGADERM+PAD 9X25 CM	CX. 25
01091000	TEGADERM+PAD 5X7 CM	CX. 50
01091001	TEGADERM+PAD 9X15 CM	CX. 25
01091002	TEGADERM+PAD 9X10 CM	CX. 25