

TEGADERM DRESSING TRANSPARENT



MAIN SPECIFICATIONS

- Tegaderm Brand transparent dressings can be used to cover and protect catheter sites and wounds, to maintain a moist environment for wound healing, and to secure devices to the skin.
- In vitro testing shows that the transparent film provides a viral barrier from viruses 27 nm in diameter or larger while the dressing remains intact without leakage.
- Waterproof, sterile barrier: protects the I.V. site or wound from external contaminants such as bacteria, viruses, blood and body fluids.
- Breathable film lets oxygen in and moisture vapor out.
- Just the right balance in adhesive strength: Tegaderm Brand dressings are made with a hypoallergenic, latex-free adhesive that is gentle to the skin yet securely holds catheters and other devices in place.
- Conformable: dressing flexes with skin for greater patient comfort.

Code	Description	Selling unit
01091010	TEGADERM DRESSING TRANSPARENT 10x12 CM	CX. 50