

ABSORBENT DRESSING STERILE NON ADHESIVE



MAIN SPECIFICATIONS

- Ideal for situations in which a high absorption capacity is needed, without adhering to the zone where they will be used.
- Because they are non-adherent, they are also used for padding areas with scabs.
- Compounds of five layers: Layer contact with the wound blank polypropylene. Layer Blue not 100% polypropylene SMS fabric with hydrophobic properties. Inner layer of cellulose pulp with high capacity absorbency. Intermediate layers on paper with spread properties. Hotmelt. Nonstick network.
- Packed and sterilized individually and then in cardboard dispenser boxes.

| Code | Description | Selling unit |
|----------|-------------------------------------|--------------|
| 01016003 | ABSORBENT DRESSING STERILE 20X20 CM | CX. 30 |
| 01016004 | ABSORBENT DRESSING STERILE 10X10 CM | CX. 60 |
| 01016018 | ABSORBENT DRESSING STERILE 10X20 CM | CX. 30 |
| 01016020 | ABSORBENT DRESSING STERILE 20X30 CM | CX. 30 |
| 01016027 | ABSORBENT DRESSING STERILE 15X20 CM | CX. 10 |

FACTORMED
Produtos Médicos & Hospitalares