

## FOAM DRESSING



### MAIN SPECIFICATIONS

- Foam covered with a film breathable and water resistant barrier and bacteria, both polyurethane materials.
- For wounds with moderate / abundant exudate.
- Highly absorbent.
- Soft and comfortable.
- Esterile.

Code	Description	Selling unit
01059001	FOAM DRESSING 10X10 CM	CX. 10
01059002	FOAM DRESSING 15X15 CM	CX. 5
01059003	FOAM DRESSING 15X20 CM	CX. 5
01059004	FOAM DRESSING 18X22 CM	CX. 5