

# ADAPTIC NON ADHERING DRESSING



## MAIN SPECIFICATIONS

- Adaptic Non-Adhering Dressing protects the fragile tissue in wounds by its unique structure with small mesh size, preventing tissue adherence to either the Adaptic dressing or the secondary dressing on top of it
- The mesh allows for exudate to easily pass through to the secondary dressing used, preventing maceration of the wound surface
- Knitted cellulose acetate fabric allows Adaptic to be cut to wound size without unravelling or shredding
- Adaptic Non-Adhering Dressing removes cleanly from the wound, minimising trauma and pain
- Adaptic Non-Adhering Dressing is a primary dressing made of knitted cellulose acetate fabric and impregnated with a specially formulated petrolatum emulsion.  
It is designed to help protect the wound while preventing the dressing from adhering to the wound<sup>1</sup> and to minimize pain and trauma upon removal.
- Adaptic Non-Adhering Dressing protects the fragile tissue in wounds by its unique structure with small mesh size, preventing tissue adherence to either the Adaptic dressing or the secondary dressing on top of it.
- The mesh allows for exudate to easily pass through to the secondary dressing used, preventing maceration of the wound surface
- Knitted cellulose acetate fabric allows Adaptic to be cut to wound size without unravelling or shredding.

- Adaptic Non-Adhering Dressing removes cleanly from the wound, minimizing trauma and pain.
- This dressing is indicated for dry to highly exuding wounds where adherence of dressing and exudate is to be prevented, including: First and second degree burns, abrasions, grafts, venous ulcers, pressure ulcers, nail extractions, eczema, staples, surgical incisions, lacerations, reconstructive procedures and suture lines.
- Adaptic non-Adhering Dressing construction Petrolatum emulsion coated mesh, Smooth and gentle, Smaller mesh size – reduced potential of secondary dressing adherence to wound bed through the mesh; Knit structure cellulose acetate.
- Can be cut to wound size.
- No loose fragments or shed on cutting.

Code	Description	Selling unit
01108007	ADAPTIC 7,6 x 7,6 CM	CX. 50
01108008	ADAPTIC 7,6 x 20,3 CM	CX. 24
01108009	ADAPTIC 7,6 x 40,6 CM	CX. 36