

## ZETUVIT E ABSORBENT DRESSING



### MAIN SPECIFICATIONS

- Hypoallergenic fibres permeable to air and water - Prevents irritation of surrounding skin and enables skin's natural functions to continue.
- Filling of fluff pulp and fluid dispersing tissue - Provide high level of absorbency.
- Blue, hydrophobic backing - Provides a barrier against contamination.
- The highly absorbent dressing with soft non-woven covering.
- Zetuvit E is a highly absorbent wound dressing pad with non-adherent, non-woven cover and absorbent cellulose fluff pad. Zetuvit E was especially developed to treat heavily exuding wounds.
- The dressing consists of 4 layers of different materials that complement one another perfectly: non-adherent non-woven cover, absorbent core of soft cellulose fluff, tissue layer for distribution of exudate and moisture-repellent backing of polypropylene non-woven used as protection against contamination.
- Available in sterile and non sterile presentation.

Code	Description	Selling unit
01046022	ZETUVIT E NON STERILE DRESSING 10X20 CM	CX. 50