

URINARY CATHETER



MAIN SPECIFICATIONS

- Soft, frosted and kink resistant PVC tubing.
- Proximal end is fitted with Universal Funnel shaped connector for safe connection to Urine bag.
- Catheters are manufactured from non-toxic, non-irritant, soft PVC compatible with Catheter lubricants.
- Available with the option of Radio opaque line throughout the length for X-ray visualisation.
- Atraumatic, soft, rounded, closed tip with two lateral eyes for efficient drainage.
- Colour coded plain connector for easy identification of size.
- Sterile product for single use.
- Length 40 cm for males and 20 cm for females.
- Available in different sizes.

Code	Description	Selling unit
09011064	URINARY MALE CATHETER 40 CM CH 10	CX.100
09011038	URINARY MALE CATHETER 40 CM CH 14	CX.100
09011040	URINARY MALE CATHETER 40 CM CH 8	CX.100
09011089	URINARY MALE CATHETER 40 CM CH 16	CX.100
09011050	URINARY FEMALE CATHETER 20 CM CH 8	CX.100

09011001	URINARY FEMALE CATHETER 20 CM CH 14	CX.100
09011065	URINARY FEMALE CATHETER 20 CM CH 10	CX.100
09011066	URINARY FEMALE CATHETER 20 CM CH 12	CX.100