

SUCTION CATHETER



MAIN SPECIFICATIONS

- Suction catheter with finger tip vacuum control.
- Non Toxic.
- Sterilised by Ethylene-oxide gas.
- Disposable device, single use.
- Frosted surface for easy insertion.
- Appropriate degree of hardness provides kink resistance.
- Atraumatic, soft rounded, bevelled open tip.
- Two lateral eyes with smooth edges.
- Color coded.

Code	Description	Selling unit
09085001	SUCTION CATHETER VACUUM CONTROL (BLACK) CH 10	CX. 100
09085006	SUCTION CATHETER VACUUM CONTROL (LIGHT GREEN) CH 06	CX. 100
09085011	SUCTION CATHETER VACUUM CONTROL (WHITE) CH 12	CX. 100

09085013	SUCTION CATHETER VACUUM CONTROL (BLUE) CH 08	CX. 100
09085014	SUCTION CATHETER VACUUM CONTROL (RED) CH 18	CX. 100
09085015	SUCTION CATHETER VACUUM CONTROL (GREEN) CH 14	CX. 100
09085016	SUCTION CATHETER VACUUM CONTROL (ORANGE) CH 16	CX. 100