

## STOMACH TUBE PVC



### MAIN SPECIFICATIONS

- Clear, medical grade PVC tube creates ease of visualization of fluids leaving and entering the stomach.
- Thermal sensitive PVC softens when warmed to body temperature for patient comfort.
- X-ray opaque line that allows easy tube and tip identification via X-ray.
- Open distal end with 4 side holes Allows for smooth suction and less tube clogging.
- Integral funnel connector is ready to attach to suction tubing.
- Available in diferent sizes.

Code	Description	Selling unit
09016000	STOMACH TUBE TYPE LEVIN RX OPAQUE CH 16 121 CM	CX. 100
09016001	STOMACH TUBE TYPE LEVIN RX OPAQUE CH 18 121 CM	CX. 100
09016006	STOMACH TUBE TYPE LEVIN RX OPAQUE CH 20 121 CM	CX.100
09016007	STOMACH TUBE TYPE LEVIN RX OPAQUE CH 14 121 CM	CX.100