

MYRINGOTOMY INSTRUMENT SET



MAIN SPECIFICATIONS

- A ready made Kit for Paracentesis Surgery.
- Single use.
- Sterile, single packed in peel-pack.
- Contents: Metal Suction Tip 1,4 x 70 mm, 2,0 x 70 mm, 2,5 x 80 mm Myringotomy Scalpel; Suction Tip Handle; Aural Speculas (2,5 mm, 4 mm, 5 mm); 2 Non woven gauzes; Sterile field.

| Code | Description | Selling unit |
|----------|---|--------------|
| 09155004 | MYRINGOTOMY BACTERIOLOGICAL EXAMINATION KIT | CX. 30 |