

# ENDOTRACHEAL TUBE REINFORCED WITH CUFF



## MAIN SPECIFICATIONS

- Designed for high volume, low pressure.
- Reinforced tube.
- Kink resistant ensures no kinking while softening at body temperature.
- Double line depth guide facilitates accurate cuff placement below vocal cord.
- Ultra-smooth material permits gentle conformity to tracheal contours.

Code	Description	Selling unit
09158003	ENDOTRACHEAL TUBES REINFORCED WITH CUFF NR 7,5	CX. 10
09158004	ENDOTRACHEAL TUBES REINFORCED WITH CUFF NR 8	CX. 10
09158002	ENDOTRACHEAL TUBES REINFORCED WITH CUFF NR 7	CX. 10