

# ABSORBING DRESSING NON STERILE NON ADHESIVE



## MAIN SPECIFICATIONS

- Ideal for situations in which a high absorption capacity is needed, without adhering to the zone where they will be used.
- Because they are non-adherent, they are also used for padding areas with scabs.
- Compounds of five layers: Layer contact with the wound blank polypropylene. Layer Blue not 100% polypropylene SMS fabric with hydrophobic properties. Inner layer of cellulose pulp with high capacity absorbency. Intermediate layers on paper with spread properties. Hotmelt. Nonstick network.
- Packed in plastic bags with 25 units.

Code	Description	Selling unit
01016031	ABSORBING DRESSING NON STERILE NON ADHESIVE 10X10 CM	CX. 25
01016032	ABSORBING DRESSING NON STERILE NON ADHESIVE 10X20 CM	CX. 25
01016033	ABSORBING DRESSING NON STERILE NON ADHESIVE 20X20 CM	CX. 25

01016034

ABSORBING DRESSING NON STERILE NON ADHESIVE 20X40  
CM

CX. 25

---