

## INTUBATION STYLET



### MAIN SPECIFICATIONS

- An aluminum compound coated rod with a plastic sheath (LDPE - low density polyethylene).
- Used to facilitate oral intubation of the trachea.
- Keep in a cool dry place.
- Individually packaged.
- Sterilized in ethylene oxide.

Code	Description	Selling unit
09158025	INTUBATION STYLET NR 5	CX. 20
09158024	INTUBATION STYLET	CX. 20
09158005	INTUBATION STYLET	CX. 20
09158026	INTUBATION STYLET	CX. 20