

OMNIFIX NON WOVEN ADHESIVE



MAIN SPECIFICATIONS

- Self-adhesive band made of a flexible soft nonwoven fabric, nonwoven fabric coated with a synthetic adhesive-free non-sensitizing synthetic latex.
- Omnifix is ideal for patients with hypersensitive skin and a tendency to allergies.
- Soft and smooth, it is easy to apply, even in the most difficult areas of application.
- Due to its high elasticity, it adapts well to body contours and will not restrict movement.
- Individually packaged.
- Soft material, non-woven.
- Synthetic adhesive secure attachment.
- Permeable to air and water vapor, however the adhesive is not sensitive to water or secretions, thus preserving its initial properties fastening, even when wet.
- It does not absorb X-rays.
- Durable and resistant to temperature changes.
- It can be removed painlessly and without leaving marks.
- Slotted support in wave paper for easy application.
- Used for fixing and retention equipment, probes, cannulas and, wound dressings (covering them completely, especially in angular joints and body parts), catheters and drainage lines.

Code	Description	Selling unit
03046024	OMNIFIX E ADHESIVE 10 M X 5 CM	CX. 1
03046019	OMNIFIX E ADHESIVE 10 M X 20 CM	CX. 1
03046004	OMNIFIX E ADHESIVE 10 M X 10 CM	CX. 1
03046005	OMNIFIX E ADHESIVE 10 M X 15 CM	CX. 1