

TEGADERM ROLL TRANSPARENT FILM



MAIN SPECIFICATIONS

- Tegaderm is a transparent roll film, not sterile, in very easy to handle roll, even with gloves.
- With two protections that promote easy application and perfect setting, with great viewing.
- Protect skin from friction, shearing and effluents.
- Will fixation waterproof, primary dressings such as gauze, alginates, foams, etc.

Code	Description	Selling unit
03091010	TEGADERM ROLL 10 M X 10 CM	CX. 1
03091011	TEGADERM ROLL TRANSPARENT 10 M X 5 CM	CX. 1