

TRANSPORE PLASTIC ADHESIVE



MAIN SPECIFICATIONS

- Easy, straight, bi-directional tear.
- Easy to handle with gloves.
- Good adhesion to skin and tubing.
- Transparent, porous.
- Hypoallergenic and not made with natural rubber latex.
- Water resistant.
- Allows clinician to tailor the tape for the various dressings and devices.
- Transparent, perforated plastic tape that tears easily, works reliably.
- Used for: securing tubing and devices (e.g., catheters, I.V. tubing), anchoring bulky dressings, when tape width must be customized.

| Code | Description | Selling unit |
|----------|------------------------------------|--------------|
| 03091004 | TRANSPORE ADHESIVE 9,14 M X1,25 CM | CX. 24 |
| 03091016 | TRANSPORE ADHESIVE 9,14 M X 2,5 CM | CX. 12 |

FACTORMED
Produtos Médicos & Hospitalares