

# HYPAFIX NON WOVEN ADHESIVE



## MAIN SPECIFICATIONS

- Hypafix consists of a white nonwoven polyester substrate coated with a skin-friendly polyacrylate adhesive, providing secure and reliable dressing retention.
- The product features a grid pattern printed on the release paper as cutting guide. The integrated split liner in the release paper eases the application of Hypafix. The soft and stretchy material allows the patient greater freedom of movement.
- Main uses: for the retention of large wound dressings, gauzes and absorbent compresses, especially on frequently moved and highly contoured parts of the body; for firm retention of instruments, probes and catheters.

Code	Description	Selling unit
03003031	HYPAFIX NON WOVEN ADHESIVE 10 M X 5 CM	CX. 1
03254017	HYPAFIX NON WOVEN ADHESIVE 10 M X 10 CM	CX. 1
03254018	HYPAFIX NON WOVEN ADHESIVE 10 M X 20 CM	CX. 1
03254002	HYPAFIX NON WOVEN ADHESIVE 10 M X 15 CM	CX. 1