

# FIXOMULL TRANSPARENT FILM



## MAIN SPECIFICATIONS

- Fixomull transparent is a waterproof wide-area fixation to reliably secure and protect wound pads and ointment dressings as well as instruments, probes, catheters and measuring devices. It is non-sterile and comes on a roll in various widths.
- Fixomull transparent consists of a thin breathable polyurethane film coated with a skin-friendly adhesive in a special pattern technique. The product protects reliably against water, bacteria and contaminants, yet granting a high MVTR. The unique application technique with the distinctive red stripe allows for easy and quick application.
- Main uses: for waterproof fixation of wound dressings, gauzes and absorbent compresses; for waterproof fixation of instruments, probes and catheters; for ointment dressings; in hydrotherapy.

| Code     | Description                      | Selling unit |
|----------|----------------------------------|--------------|
| 03254001 | FIXOMULL TRANSPARENT 10 M X15 CM | CX. 1        |
| 03254000 | FIXOMULL TRANSPARENT 10 M X10 CM | CX. 1        |
| 3254039  | FIXOMULL TRANSPARENT 10 CM X 2 M | CX. 1        |
| 72216-00 | FIXOMULL TRANSPARENT 5 CM X 10 M | CX. 1        |