

# LEUKOFLEX ADHESIVE



## MAIN SPECIFICATIONS

- Leukoflex is a transparent surgical tape ideal for securing dressings on patients with sensitive skin. It is impermeable to water, allowing patients with wound dressings secured with Leukoflex to shower, swim or participate in hydrotherapy.
- Suitable for use as an occlusive dressing to seal in topical medication, Leukoflex features a hypoallergenic adhesive which is skin-friendly and provides excellent initial and long-term adhesion. Removal of the tape is painless and does not leave residue.
- Constructed from a very thin and stretchable polyethylene film, Leukoflex conforms smoothly to body contours and adapts well to patient movement.
- Waterproof tape for impermeable dressings.
- Leukoflex is made of a transparent polyethylene film coated with a skin-friendly polyacrylate adhesive. The thin film is impermeable to moisture allowing the patients to wash or shower. The stretchable film conforms smoothly to body contours and adapts well to movements.
- Main uses: Fixation of wound dressings of all types and sizes; Secure fixation also under water For occlusive dressings.

Code	Description	Selling unit
03254031	LEUKOFLEX ADHESIVE 5 M X 5 CM	CX. 6
03254028	ADESIVO LEUKOFLEX 5X2,5 CM	CX. 12
03254038	ADESIVO LEUKOFLEX 5X1,25 CM	CX. 24