

BEQUILLE CATHETER



MAIN SPECIFICATIONS

- Bequille catheter has 35 cm in length, varying diameters.
- Set composed of catheters 16 of different sizes in Charrier scale where one unit corresponds to 1/3 mm.
- The one way catheter has a self-restraint system requiring external fixation.
- The catheter tip of angled Coudé (bequille) are useful in cases of stenosis and tightening of the prostatic urethra. Used for catheterization of the bladder evacuator in patients with prostate hyperplasia.

Code	Description	Selling unit
11158008	BEQUILLE CATHETER CH 18	CX. 5
11158005	BEQUILLE CATHETER CH 12	CX. 5
11158007	BEQUILLE CATHETER CH 16	CX. 5
11158004	BEQUILLE CATHETER CH 10	CX. 5
11158002	BEQUILLE CATHETER CH 8	CX. 5