

# FOLEY CATHETER SILICONE COATED 2 WAY



## MAIN SPECIFICATIONS

- Algalia Foley composed of latex with a silicone elastomer-based coating made from a hydrophobic material that rejects water.
- The silicone elastomer coating offers protection against urethral irritation.
- Innovative manufacturing process provides the maximum amount of silicone coating available.
- Latex catheter internally and externally coated in 100% silicone to withstand fouling.
- Symmetrical balloon for greater reliability and good positioning in the bladder.

Code	Description	Selling unit
11011001	FOLEY CATHETER SILICONE COATED 2 WAYS CH 10	CX. 10
11254000	FOLEY CATHETER SILICONE COATED 2 WAYS CH 12	CX. 10
11254001	FOLEY CATHETER SILICONE COATED 2 WAYS CH 14	CX. 10
11254002	FOLEY CATHETER SILICONE COATED 2 WAYS CH 16	CX. 10
11254003	FOLEY CATHETER SILICONE COATED 2 WAYS CH 18	CX. 10
11254005	FOLEY CATHETER SILICONE COATED 2 WAYS CH 20	CX. 10

11254006	FOLEY CATHETER SILICONE COATED 2 WAYS CH 22	CX. 10
11011008	FOLEY CATHETER SILICONE COATED 2 WAYS CH 24	CX. 10
11011009	FOLEY CATHETER SILICONE COATED 2 WAYS CH 26	CX. 10
11011014	FOLEY CATHETER SILICONE COATED 2 WAYS CH 08	CX. 10
11011023	FOLEY CATHETER SILICONE COATED 2 WAYS CH 06	CX. 10