

FOLEY CATHETER 100% SILICONE 2 WAY



MAIN SPECIFICATIONS

- Is a 100% silicone 2-way biocompatible indwelling foley catheter designed to eliminate the risk of allergic reactions and provide maximum long-term patient comfort.
- Smooth Rounded Tip - Allows for easy insertion and provides maximum patient compliance.
- Large Twin Catheter Eyes - The smooth moulded drainage eyes reduce the risk of snagging, ensuring easier insertion when in place. The large eyes allow quick fluid drainage and reduce risk of blockages.
- Large Inner Lumen - has a large inner lumen with a round, smooth surface reducing the risk of bacterial colonization and encrustation and allowing quick drainage of fluids and any possible clots.
- Recessed Inflation Valve - Prevents accidental emptying of the catheter balloon.
- Spherical Retention Balloon - The ribbed design of the balloon gives perfect symmetry providing a good position and seal when in situ. The silicone balloon has the ability to contract back to its unfilled streamline shape providing easier removal even after 30 days of use.
- Flexible Shaft - Designed to provide optimum stiffness for easy insertion whilst maintaining true flexibility, giving users freedom of movement without discomfort.
- 100% Silicone - is made of 100% silicone which minimizes bacterial colonization and encrustation whilst in situ, facilitating easier removal.
- Ribbed Outlet Connector - Provides a secure attachment to a catheter valve or leg bag, minimizing the risk of slippage, disconnection or leakage.
- Clearly marked with sizing information.

Code	Description	Selling unit
11085000	FOLEY CATHETER SILICONE 2 WAYS CH 12	CX. 10
11279000	FOLEY CATHETER SILICONE 2 WAYS CH 14	CX. 10
11279001	FOLEY CATHETER SILICONE 2 WAYS CH 16	CX. 10
11279002	FOLEY CATHETER SILICONE 2 WAYS CH 18	CX. 10
11279003	FOLEY CATHETER SILICONE 2 WAYS CH 20	CX. 10
11085005	FOLEY CATHETER SILICONE 2 WAYS CH 22	CX. 10
11085006	FOLEY CATHETER SILICONE 2 WAYS CH 24	CX. 10