

FOLEY CATHETER SILICONE COATED 3 WAY



MAIN SPECIFICATIONS

- Algalia Foley composed of latex with a silicone elastomer-based coating made from a hydrophobic material that rejects water.
- The silicone elastomer coating offers protection against urethral irritation.
- Innovative manufacturing process provides the maximum amount of silicone coating available.
- Latex catheter internally and externally coated in 100% silicone to withstand fouling.
- Symmetrical balloon for greater reliability and good positioning in the bladder.

Code	Description	Selling unit
11158011	FOLEY CATHETER SILICONE COATED 3 WAYS CH 16	CX. 10
11011022	FOLEY CATHETER SILICONE COATED 3 WAYS CH 14	CX. 10
11158012	FOLEY CATHETER SILICONE COATED 3 WAYS CH 18	CX. 10