

## URINE BAG STERILE WITH TAP



### MAIN SPECIFICATIONS

- Drainable urine bags provide a robust and versatile urine drainage solution for both hospital and community environments. All urine drainable bags are fitted with a non-return valve as standard to prevent urine reflux and protect the patient against ascending infection.
- Easy to Operate Leak-Free Taps - Both the lever tap and T-tap fitted to the bags are 100% leak tested and are easy to operate. The lever tap requires a lower force to open and close the tap and this design is often preferred by the elderly or patients with impaired dexterity.
- Free and Rapid Flow - The wide bore drainage tube allows unimpaired flow of urine into the bag, optimising drainage and providing observation of urine flow. The tube is both firm and flexible affording comfort to the user.
- Needle-Free Sample Port - The needle-free sample port eliminates the risk of needle stick injury, whilst ensuring the bag is leak free.
- Available in different sizes.

Code	Description	Selling unit
11011021	URINE BAG 1,5 L STERILIZED WITH NON RETURN VALVE AND 90 CM TUBE	CX.100
11011011	URINE BAG 1,5 L STERILIZED WITH NON RETURN VALVE, TAP AND 90 CM TUBE	CX.100