

SILICON URINARY MALE EXTERNAL CATHETER



MAIN SPECIFICATIONS

- The collector / external urinary catheter (length up to 7.5 cm) is available in five diameters to provide greater comfort to the user, better fit and greater confidence (diameter 25 mm to 41 mm).
- Latex free.
- Made of silicon: hypoallergenic, with less risk of maceration and irritation due to the ease of gas exchange.
- The adhesive strip provides easy application.
- Transparent, allows skin monitoring.
- It allows the practice of urinary catheterization.
- Compatible with all collectors bags.
- Risks of reduced infection: the torque tube allows for quick output of urine preventing building up the level of the urinary meatus.

Code	Description	Selling unit
11001003	SILICON URINARY MALE EXTERNAL CATHETER 36 MM	CX. 1

