

# PLASTIC CATHETER CLAMP



## MAIN SPECIFICATIONS

- The Clamp is a molded plastic catheter tube clamp designed to keep your catheter securely blocked to prevent leakage.
- Made with a strong, durable plastic material, with a positive lock mechanism and is easily operated one-handedly.
- Plastic catheter clamps are designed for security and durability.
- The easily-operated jaws are serrated and feature a positive lock clamp for a secure hold.
- Also, they are conveniently easy to use with one hand.
- Applications: Keeps catheters securely blocked to prevent leakage during bag changes.
- Catheter Clamp Features: Strong, durable plastic material, Serrated jaws, Positive lock clamp One-handed operation.
- Single use.
- Packed individually and sterilized by ethylene oxide.

Code	Description	Selling unit
11011018	UMBILICAL COLD PLASTIC CLAMP BLUE FOR CATHETERS	CX.100