

OXYGEN CATHETER



MAIN SPECIFICATIONS

- Purpose: Oxygen administration, replaces the use of masks.
- Available in two types: open atraumatic tip with rounded edges, with side holes 8 and foam pad.
- Medical grade PVC tube, non toxic and non-pyrogenic material, with very smooth surface that minimizes friction.
- Wireless RX, conical connector.
- Individualized packaging.

Code	Description	Selling unit
09011026	OXYGEN CATHETER WITH DETACHABLE FOAM 40 CM CH 08 (BLUE)	CX.100
09011027	OXYGEN CATHETER WITH DETACHABLE FOAM 40 CM CH 12 (WHITE)	CX.100
09011041	OXYGEN CATHETER WITH DETACHABLE FOAM 40 CM CH 14 (GREEN)	CX.100
09011063	OXYGEN CATHETER WITH DETACHABLE FOAM 40 CM CH 10	CX.100