

SPIGOT FOR CATHETERS



MAIN SPECIFICATIONS

- The catheter plug can be used to securely close the urine catheter.
- Catheter Spigot is used to provide stop of flow for catheters hygienically during nursing procedures. It is non-invasive has been utilized to seal the catheter for a short time to allow urine to collect in the bladder.
- Universal size is molded to fit all sizes of catheters or tubes (conical shape, adapts to any size or different size tubing).
- Graduated design made it to easily accommodate tubing of various internal diameters.
- Ribbed designed provides strong, secure connections.
- Top with a retaining loop prevents it being lost whilst not in use.
- Ergonomic flanged designed assists withdrawal and maximize grip.
- Can be supplied sterile in individual peel pouches, or even in bulk packaging.

Code	Description	Selling unit
11011012	SPIGOT FOR CATHETERS WITHOUT HANDLE NON STERILE	CX. 100
11011015	SPIGOT FOR CATHETERS WITH HANDLE NON STERILE	CX. 100

FACTORMED
Produtos Médicos & Hospitalares