

# LATEX URINARY MALE EXTERNAL CATHETER



## MAIN SPECIFICATIONS

- Penile Sheath/External Catheter.
- Is specially designed for urine incontinence for day and night use in male patients.
- Manufactured from 100% latex for soft and gentle feel.
- Provided with soft adhesive coated strip for proper fixing on the penis.
- Proximal end is designed for safe and easy connection to urine bag/leg bag.
- Unwind easily, facilitating its application.

Code	Description	Selling unit
11034000	URINARY MALE EXTERNAL CATHETER MEDIUM 25 MM	CX. 50
11034001	URINARY MALE EXTERNAL CATHETER LARGE 30 MM	CX. 50
11034002	URINARY MALE EXTERNAL CATHETER XL 35 MM	CX. 50