

## CATGUT SUTURE



### MAIN SPECIFICATIONS

- High tensile strength.
- Excellent tissue support.
- Absorption by phagocytosis.
- Uniform absorption.
- Twisted structure.
- Excellent handling properties.
- High flexibility.
- High tensile strength.
- Excellent knot security.
- Uniform diameter by centerless polishing (monofilament appearance).
- Soft passage through the tissues without sawing, tissue drag and trauma.
- Excellent knot tie down.
- Hermetically sealed packing.
- Guaranteed seal and product sterility.
- Indications: General soft tissue closing and/or ligation especially in General surgery; Gastrointestinal

surgery; Gynaecology; Obstetrics; Urology; Ophthalmic surgery.

- Contraindications: Cardiovascular surgery Neurosurgery. Being absorbable not to use when an extended wound support is required. Caution in patients with delayed wound healing (e.g. compromised immune system). In the region of infected wounds Known sensitivities or allergies to collagen.

Code	Description	Selling unit
08113000	CATGUT CHROM 1 75 CM DS37 LANCEOLATE NEEDLE 37 MM	CX. 12
08113001	CATGUT CHROM 2/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113002	CATGUT CHROM 0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113024	CATGUT CHROM 2/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113020	CATGUT CHROM 5/0 45 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113081	CATGUT PLAIN 1 75 CM HR37 CYLINDER NEEDLE 37 MM	CX. 12
08113082	CATGUT PLAIN 2 75 CM HR45 CYLINDER NEEDLE 45 MM	CX. 12
08113083	CATGUT CHROM 1 75 CM HR48 CYLINDER NEEDLE 48 MM	CX. 12
08113084	CATGUT CHROM 2 75 CM HR45 CYLINDER NEEDLE 45 MM	CX. 12
08113032	CATGUT CHROM 0 75 CM DS37 AG/LANCEOLADA 37 MM	CX. 12
08113213	CATGUT CHROM 1 75 CM DS24 AG/LANCEOLADA 24 MM	CX. 12
08113215	CATGUT CHROM 3/0 75 CM DS24 AG/LANCEOLADA 24 MM	CX. 12
08113216	CATGUT CHROM 4/0 75 CM DS24 AG/LANCEOLADA 24 MM	CX. 12