

SURGICRYL PGA SUTURE POLYGLYCOLID ACID



MAIN SPECIFICATIONS

- High tensile strength.
- Excellent tissue support.
- Synthetic Minimal tissue reaction.
- Breaks down by hydrolysis.
- Predictable absorption rate
- Braided structure.
- Excellent handling properties.
- High flexibility.
- High tensile strength.
- Excellent knot security.
- Coated multifilament.
- Soft passage through the tissues with minimal sawing, tissue drag and trauma.
- Good knot tie down / adjustability.
-

Reduced capillary action.

- Hermetically sealed packing.
- Guaranteed seal and product sterility.
- Indications: General soft tissue closing and/or ligation especially in: General surgery, Skin closure, Gastrointestinal surgery, Gynaecology Obstetrics, Plastic surgery, Urology, Ophthalmic surgery, Orthopaedics.
- Contraindications: Cardiovascular surgery, Neurosurgery. Being absorbable not to use when an extended wound support is required. Caution in patients with delayed wound healing (e.g. compromised immune system).

Code	Description	Selling unit
08113003	SURGICRYL PGA VIOLET 3/0 75 CM HR22 CYLINDER NEEDLE 22 MM	CX. 12
08113004	SURGICRYL PGA VIOLET 4/0 75 CM HR22 CYLINDER NEEDLE 22 MM	CX. 12
08113005	SURGICRYL PGA VIOLET 4/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113006	SURGICRYL PGA VIOLET 3/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113007	SURGICRYL PGA VIOLET 2/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113011	SURGICRYL PGA VIOLET 0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113012	SURGICRYL PGA VIOLET 0 75 CM DS30 LANCEOLATE NEEDLE 30 MM	CX. 12
08113013	SURGICRYL PGA VIOLET 0 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113014	SURGICRYL PGA VIOLET 1 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113015	SURGICRYL PGA VIOLET 2/0 75 CM DS30 CYLINDER NEEDLE 30 MM	CX. 12
08113016	SURGICRYL PGA VIOLET 5/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113017	SURGICRYL PGA VIOLET 1 75 CM DS36 LANCEOLATE NEEDLE 36 MM	CX. 12
08113021	SURGICRYL PGA VIOLET 1 75 CM DS40 LANCEOLATE NEEDLE 40 MM	CX. 12

08113027	SURGICRYL PGA VIOLET 1 75 CM HR36 CYLINDER NEEDLE 36 MM	CX. 12
08113028	SURGICRYL PGA VIOLET 6/0 45 CM DS12 LANCEOLATE NEEDLE 12 MM	CX.12
08113034	SURGICRYL PGA VIOLET 0 90 CM HR36 CYLINDER NEEDLE 36 MM	CX. 12
08113035	SURGICRYL PGA VIOLET 0 75 CM DS36 LANCEOLATE NEEDLE 36 MM	CX. 12
08113043	SURGICRYL PGA VIOLET 2/0 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113050	SURGICRYL PGA VIOLET 3/0 75 CM HR17 CYLINDER NEEDLE 17 MM	CX. 12
08113051	SURGICRYL PGA VIOLET 2/0 75 CM HR22 CYLINDER NEEDLE 22 MM	CX. 12
08113058	SURGICRYL PGA VIOLET 3/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113059	SURGICRYL 910 VIOLET 4/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113139	SURGICRYL PGA VIOLET 3/0 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113155	SURGICRYL 910 VIOLET 3/0 75 CM HR30 CYLINDER NEEDLE 30 MM	CX. 12
08113162	SURGICRYL PGA VIOLET 2 90 CM HR48 CYLINDER NEEDLE 48 MM	CX. 12
08113168	SURGICRYL PGA VIOLET 3/0 75 CM HS26 LANCEOLATE NEEDLE 26MM	CX. 12
08113169	SURGICRYL PGA VIOLET 5/0 75 CM DS16 LANCEOLATE NEEDLE 16 MM	CX. 12
08113177	SURGICRYL PGA VIOLET 2/0 75 CM HR30 CYLINDER NEEDLE 30 MM	CX. 12
08113093	SURGICRYL PGA VIOLET 1 90 CM HS36 LANCEOLATE NEEDLE 36 MM	CX. 12
08113152	SURGICRYL PGA VIOLET 2/0 75 CM HR36 CYLINDER NEEDLE 36 MM	CX. 12
08113153	SURGICRYL PGA VIOLET 0 75 CM HR30 CYLINDER NEEDLE 30 MM	CX. 12
08113176	SURGICRYL PGA VIOLET 1 90 CM HS40 LANCEOLATE NEEDLE 40 MM	CX. 12
08113197	SURGICRYL PGA VIOLETA 0 90 CM HR48 AG/LANCEOLADA 1/2 48 MM	CX. 12

08113196	SURGICRYL PGA VIOLETA 0 90 CM HS48 AG/LANCEOLADA 1/2 48 MM	CX. 12
8113093	SURGICRYL PGA VIOLETA 1 90 CM HS36 AG/LANCEOLADA 36 MM	CX. 12
08113199	SURGICRYL PGA VIOLETA 3/0 75 CM DR17 AG/LANCEOLADA 17 MM CX.12	CX. 12
08113217	SURGICRYL PGA VIOLETA 4/0 75 CM DS16 AG/LANCEOLADA 16 MM	CX. 12
08113193	SURGICRYL PGA VIOLETA 4/0 75 CM HR17 AG/CILINDRICA 17 MM	CX. 12
08113200	SURGICRYL PGA VIOLETA 7/0 45 CM 2 X 8 mm - 200µm ½ CIRCULO	CX. 12
08113096	CASSETE SURGICRYL PGA VIOLETA 3/0 E 4/0 15 MM	CX. 12
08113105	CASSETE SURGICRYL PGA VIOLETA 6/0 15 MM	UN