

SILK SUTURE



MAIN SPECIFICATIONS

- High tensile strength, non absorbable
- Good and extended tissue support for up to 3 month.
- Braided or twisted structure
- Excellent handling properties.
- High flexibility.
- Excellent knot security.
- Coated multifilament.
- Soft passage through the tissues with minimal sawing, tissue drag and trauma.
- Good knot tie down / adjustability.
- Reduced capillary action.
- Hermetically sealed packing
- Guaranteed seal and product sterility.
- Indications: General soft tissue closing and/or ligation especially in: General surgery, Skin closure, Ophthalmic surgery, Oral surgery, Gastrointestinal surgery, Gynaecology, Obstetrics.
- Contraindications: Urinary tract, Biliary tract, Known allergies or sensitivities to silk.

Code	Description	Selling unit
08113008	SILK BRAIDED BLACK 2/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113009	SILK BRAIDED BLACK 3/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113018	SILK BRAIDED BLACK 4/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113019	SILK BRAIDED BLACK 5/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113022	SILK BRAIDED BLACK 4/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113029	SILK BRAIDED BLACK 0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113030	SILK BRAIDED BLACK 1 75 CM DS30 LANCEOLATE NEEDLE 30 MM	CX. 12
08113031	SILK BRAIDED BLACK 0 75 CM DS30 LANCEOLATE NEEDLE 30 MM	CX. 12
08113037	SILK BRAIDED BLACK 3/0 75 CM HR20 CYLINDER NEEDLE 20 MM	CX. 12
08113038	SILK BRAIDED BLACK 2/0 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113044	SILK VIRGIN BLACK 6/0 45 CM 2 X 6,5 MM - 240 MCM	CX. 12
08113045	SILK VIRGIN BLACK 6/0 45 CM 1/4 2 X 8 MM - 250 MCM	CX. 12
08113048	SILK BRAIDED BLACK 5/0 45 CM DS12 LANCEOLATE NEEDLE 12 MM	CX. 12
08113056	SILK BRAIDED BLACK 5/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113057	SILK BRAIDED BLACK 3/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113068	SILK BRAIDED BLACK 0 75 CM HS36 LANCEOLATE NEEDLE 36 MM	CX. 12
08113069	SILK BRAIDED BLACK 5/0 75 CM DS16 LANCEOLATE NEEDLE 16 MM	CX. 12
08113070	SILK BRAIDED BLACK 6/0 45 CM DS12 LANCEOLATE NEEDLE 12 MM	CX. 12
08113076	SILK BRAIDED BLACK 3/0 75 CM DS16 LANCEOLATE NEEDLE 16 MM	CX. 12
08113077	SILK BRAIDED BLACK 4/0 75 CM HS17 LANCEOLATE NEEDLE 17	CX. 12

MM

08113078	SILK BRAIDED BLACK 2/0 75 CM HS26 LANCEOLATE NEEDLE 26 MM	CX. 12
08113100	SILK BRAIDED BLACK 6/0 75 CM DS16 LANCEOLATE NEEDLE 16 MM	CX. 12
08113107	SILK BRAIDED BLACK 2/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113112	SILK BRAIDED BLACK 2/0 75 CM HR20 CYLINDER NEEDLE 20 MM	CX. 12
08113111	SILK BRAIDED BLACK 1 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113113	SILK BRAIDED BLACK 1/0 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113114	SILK BRAIDED BLACK 2/0 75 CM DS30 LANCEOLATE NEEDLE 30 MM	CX. 12
08113115	SILK BRAIDED BLACK 2/0 75 CM DS16 LANCEOLATE NEEDLE 16 MM	CX. 12
08113116	SILK BRAIDED BLACK 3/0 75 CM DS30 LANCEOLATE NEEDLE 30 MM	CX. 12
08113117	SILK BRAIDED BLACK 3/0 75 CM HS26 LANCEOLATE NEEDLE 26 MM	CX. 12
08113118	SILK BRAIDED BLACK 4/0 75 CM DS16 LANCEOLATE NEEDLE 16 MM	CX. 12
08113119	SILK BRAIDED BLACK 4/0 75 CM HR20 CYLINDER NEEDLE 20 MM 1/2 CIRCLE	CX. 12
08113120	SILK BRAIDED BLACK 3/0 75 CM DS26 LANCEOLATE NEEDLE 26 MM 1/2 CIRCLE	CX. 12
08113121	SILK BRAIDED BLACK 3/0 75 CM HS17 LANCEOLATE NEEDLE 17 MM 1/2 CIRCLE	CX. 12
08113123	SILK BRAIDED BLACK 2 75 CM DS37 LANCEOLATE NEEDLE 37 MM 3/8	CX. 12
08113130	SILK BRAIDED BLACK 2/0 75 CM DS26 LANCEOLATE NEEDLE 26 MM	CX. 12
08113134	SILK BRAIDED BLACK 0 75 CM HR30 CYLINDER NEEDLE 30 MM	CX. 12