

SUPRAMID SUTURE POLYAMIDE



MAIN SPECIFICATIONS

- High tensile strength, non absorbable.
- Excellent and permanent tissue support.
- Synthetic.
- Minimal tissue reaction.
- Twisted structure.
- Excellent handling properties.
- High flexibility.
- Excellent knot security.
- Multifilament in a sheath.
- Soft passage through the tissues without sawing, tissue drag and trauma.
- Good knot tie down / adjustability.
- Reduced capillary action.
- Pliability and softness.
- Easy and secure knotting.
- Hermetically sealed packing.

- Guaranteed seal and product sterility.
- Indications: General soft tissue closing and/or ligation especially in skin closure.

Code	Description	Selling unit
08113023	SUPRAMID BLACK 2/0 75 CM DS24 LANCEOLATE NEDDLE 24 MM	CX. 12
08113026	SUPRAMID BLACK 2/0 75 CM DS19 LANCEOLATE NEDDLE 19 MM	CX. 12
08113039	SUPRAMID BLACK 0 75 CM DS24 LANCEOLATE NEDDLE 24 MM	CX. 12
08113161	SUPRAMID BLACK 1 75 CM DS36 LANCEOLATE NEEDLE 36 MM	CX. 12
08113067	SUPRAMID BLACK 3/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113101	SUPRAMID BLACK 4/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 36
08113092	SET SUPRAMID BLACK 4/0 100 MM	CX. 1
08113142	SET SUPRAMID BLACK 3&4 50 MM	CX. 1
08001013	SUPRAMID BLACK 5/0 75CM DS19 AG/LANCEOLADA 3/8 19MM CX.36	CX. 36
08001016	SUPRAMID BLACK 1 75CM HS26s AG/LANCEOLADA 1/2 26 MM	CX. 36