

# SURGICRYL RAPID SUTURE POLYGLYCOLID ACID MULTIFILAMENT



## MAIN SPECIFICATIONS

- High initial tensile strength.
- Excellent tissue support.
- Synthetic.
- Minimal tissue reaction.
- Breaks down by hydrolysis.
- Predictable absorption rate, quick mass absorption.
- Braided structure.
- Excellent handling properties.
- High flexibility.
- High tensile strength.
- Excellent knot security.
- Coated multifilament.
- Soft passage through the tissues with minimal sawing, tissue drag and trauma.
-

Good knot tie down / adjustability.

- Reduced capillary action.
- Hermetically sealed packing.
- Guaranteed seal and product sterility.
- Indications: Superficial soft tissue, closing of skin and mucosa only.
- especially in Plastic surgery, Obstetrics (Episiotomy repair), Gynaecology, Urology (Circumcision), Ophthalmic surgery (Conjunctiva) and Pediatric surgery.
- Contraindications: Cardiovascular surgery, Neurosurgery. When wound support for more than 7 days is required ligation. Caution in patients with delayed wound healing (e.g. compromised immune system).

Code	Description	Selling unit
08113041	SURGICRYL RAPID 5/0 45 CM DS12 LANCEOLATE NEEDLE 3/8 12 MM	CX. 12
08113064	SURGICRYL RAPID 3/0 75 CM DS24 LANCEOLATE NEEDLE 3/8 24 MM	CX. 12
08113128	SURGICRYL RAPID 3/0 150 CM SEM AGULHA WITHOUT NEEDLE	CX. 12
08113131	SURGICRYL RAPID 2/0 75 CM DS24 LANCEOLATE NEEDLE 3/8 24 MM	CX. 12
08113140	SURGICRYL RAPID 3/0 45 CM DS19 LANCEOLATE NEEDLE 3/8 19 MM	CX. 12
08113164	SURGICRYL RAPID 4/0 75 CM DS19 LANCEOLATE NEEDLE 3/8 19 MM	CX. 12
08113165	SURGICRYL RAPID 6/0 75 CM DS12 LANCEOLATE NEEDLE 3/8 12 MM CX.12	CX. 12
08113166	SURGICRYL RAPID 4/0 75 CM HR22 CYLINDER NEEDLE 1/2 22 MM	CX. 12
08113167	SURGICRYL RAPID 5/0 75 CM HR17 CYLINDER NEEDLE 1/2 17 MM	CX. 12
08113186	SURGICRYL RAPID 3/0 45 CM HS17 AG/LANCEOLADA 1/2 17 MM	CX. 12
08113203	SURGICRYL RAPID 4/0 75 CM DS16 AG/LANCEOLADA 3/8 16 MM	CX. 12
08113207	SURGICRYL RAPID 5/0 50 CM DS12 AG/LANCEOLADA 3/8 EASY PASS 12 MM	CX. 12
08113198	SURGICRYL RAPID 5/0 75 CM DS16 AG/LANCEOLADA 3/8 16 MM	CX. 12