

SURGICRYL MONOFAST SUTURE POLYGLYCAPRONE 25 MONOFILAMENT



MAIN SPECIFICATIONS

- High tensile strength.
- Excellent tissue support.
- Synthetic.
- Minimal tissue reaction.
- Breaks down by hydrolysis Predictable absorption rate.
- Monofilament with microscopic uniform diameter.
- Soft passage through the tissues without sawing, tissue drag and trauma.
- Absence of capillarity (no support of bacterial growth).
- Excellent knot tie down.
- Pliability and softness.
- Good handling properties.
- Easy and secure knotting.
- Hermetically sealed packing.
-

Guaranteed seal and product sterility.

- Indications: Superficial soft tissue and/or ligation especially in General surgery, Skin closure, Plastic surgery, Gastrointestinal surgery, Gynaecology, Urology.
- Contraindications: Cardiovascular surgery, Neurosurgery, Ophthalmic surgery. Being absorbable not to use when an extended wound support is required. Caution in patients with delayed wound healing (e.g. compromised immune system).

Code	Description	Selling unit
08113042	SURGICRYL MONOFAST COLORLESS 2/0 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113073	SURGICRYL MONOFAST COLORLESS 2/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113122	SURGICRYL MONOFAST COLORLESS 4/0 45 CM DS12 LANCEOLATE NEEDLE 12 MM	CX. 12
08113150	SURGICRYL MONOFAST COLORLESS 5/0 45 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113179	SURGICRYL MONOFAST VIOLET 2/0 75 CM HR30 CYLINDER NEEDLE 30 MM	CX. 12