

POLYESTER SUTURE



MAIN SPECIFICATIONS

- Unsurpassed tensile strength, non absorbable.
- Excellent and permanent tissue support.
- Resistance to repeated bending.
- Use in cardiovascular surgery (heart prosthesis).
- Synthetic.
- Minimal tissue reaction.
- Braided structure.
- Excellent handling properties.
- High flexibility.
- Excellent knot security.
- Coated multifilament.
- Soft passage through the tissues with minimal sawing, tissue drag and trauma.
- Good knot tie down / adjustability.
- Reduced capillary action.
- Hermetically sealed packing.

- Guaranteed seal and product sterility.
- Indications: General soft tissue closing and/or ligation especially in Cardiovascular surgery, Ophthalmic surgery, General surgery, Skin closure, Orthopaedics, Gastrointestinal surgery, Obstetrics.

Code	Description	Selling unit
08113046	POLYESTER GREEN 1 75 CM DS30 LANCEOLATE NEEDLE 30 MM	CX. 12
08113173	POLYESTER GREEN 2/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113174	POLYESTER GREEN 5/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113175	POLYESTER GREEN 3/0 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113204	POLYESTER GREEN 1 75 CM HS36 ½ AG/LANCEOLADA 36 MM	CX. 12