

SURGICRYL 910 SUTURE POLYGLACTINE



MAIN SPECIFICATIONS

- High tensile strength.
- Excellent tissue support.
- Synthetic.
- Minimal tissue reaction.
- Breaks down by hydrolysis.
- Predictable absorption rate.
- Braided structure.
- Excellent handling properties.
- High flexibility.
- Excellent knot security.
- Coated multifilament.
- Soft passage through the tissues with minimal sawing, tissue drag and trauma.
- Good knot tie down / adjustability.
- Reduced capillary action.
- Hermetically sealed packing.

- Guaranteed seal and product sterility.
- Indications: General soft tissue closing and/or ligation especially in General surgery, Skin closure, Gastrointestinal surgery, Plastic surgery, Ophthalmic surgery, Gynaecology, Obstetrics, Orthopaedics and Urology.
- Contraindications: Cardiovascular surgery, Neurosurgery. Being absorbable not to use when an extended wound support is required. Caution in patients with delayed wound healing (e.g. compromised immune system).

Code	Description	Selling unit
08113059	SURGICRYL 910 VIOLET 4/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113060	SURGICRYL 910 VIOLET 5/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113061	SURGICRYL 910 VIOLET 3/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113072	SURGICRYL 910 VIOLET 2/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113079	SURGICRYL 910 VIOLET 3/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113080	SURGICRYL 910 VIOLET 3/0 75 CM DS19 LANCEOLATE NEEDLE 19 MM	CX. 12
08113097	SURGICRYL 910 VIOLET 0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113098	SURGICRYL 910 VIOLET 3/0 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113099	SURGICRYL 910 VIOLET 2 75 CM HR36 CYLINDER NEEDLE 36 MM	CX. 12
08113127	SURGICRYL 910 VIOLET 0 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113137	SURGICRYL 910 VIOLET 2/0 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113138	SURGICRYL 910 VIOLET 6/0 75 CM DS12 CYLINDER NEEDLE 12 MM	CX. 12
08113143	SURGICRYL 910 VIOLET 2/0 75 CM HR40 CYLINDER NEEDLE 40 MM	CX. 12
08113144	SURGICRYL 910 VIOLET 2 90 CM HR40 CYLINDER NEEDLE 40 MM	CX. 12

08113154	SURGICRYL 910 VIOLET 2/0 75 CM HR30 CYLINDER NEEDLE 30 MM	CX. 12
08113155	SURGICRYL 910 VIOLET 3/0 75 CM HR30 CYLINDER NEEDLE 30 MM	CX. 12
08113202	SURGICRYL 910 VIOLETA 1 75 CM HR30 AG/CILINDRICA 30 MM	CX. 12
08113188	SURGICRYL 910 VIOLETA 4/0 75 CM HS17 AG/LANCEOLADA 1/2 17 MM	CX. 12
08113185	SURGICRYL 910 VIOLETA 7/0 45 CM 2 X 8 mm - 200µm 1/2 CIRCULO	CX. 12
08113208	SURGICRYL 910 INCOLOR 3/0 50 CM DS16 AG/LANCEOLADA 3/8 EASY PASS 16 MM	CX. 12