

SURGICRYL MONOFILAMENT SUTURE



MAIN SPECIFICATIONS

- High and long-lasting tensile strength.
- Outstanding strength retention for extended wound support.
- Synthetic, is a synthetic absorbable sterile surgical suture composed of the polyester poly (p-dioxanone).
- Minimal tissue reaction.
- Breaks down by hydrolysis.
- Predictable absorption rate.
- Monofilament with microscopic uniform diameter.
- Soft passage through the tissues without sawing, tissue drag and trauma.
- Absence of capillarity (no support of bacterial growth).
- Excellent knot tie down.
- Pliability and softness.
- Good handling properties
- Easy and secure knotting.
- Hermetically sealed packing.
-

Guaranteed seal and product sterility.

- Indications: General soft tissue closing especially in General surgery, Gastrointestinal surgery, Ophthalmic surgery, Orthopaedics, Abdominal wall closure, Gynaecology, Plastic surgery and Urology.
- Contraindications: Cardiovascular surgery, Neurosurgery, Prosthetic devices. Being absorbable not to use when wound support for more than 6 weeks is required. Caution in patients with delayed wound healing (e.g. compromised immune system).
- Is available in a range of gauge sizes and lengths, non-needed or attached to atraumatic stainless steel needles of varying types.

Code	Description	Selling unit
08113094	SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 2/0 75 CM HR26 CYLINDER NEEDLE 26 MM	CX. 12
08113095	SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 3/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113129	SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 0 90 CM HR48 CYLINDER NEEDLE 1/2 48 MM	CX. 12
08113158	SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 2/0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113159	SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 0 75 CM DS24 LANCEOLATE NEEDLE 24 MM	CX. 12
08113160	SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 1 75 CM HS36 LANCEOLATE NEEDLE 36 MM	CX. 12
08113191	SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLETA 0 90CM HR36 AG/CILINDRICA 1/2 36 MM	CX. 12
08113190	SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLETA 1/0 90CM HR36 AG/CILINDRICA 1/2 36 MM	CX. 12