

# SURGICRYL MONOFILAMENT SUTURE



## MAIN SPECIFICATIONS

- High and long-lasting tensile strength.
- Outstanding strength retention for extended wound support.
- Synthetic, is a synthetic absorbable sterile surgical suture composed of the polyester poly (p-dioxanone).
- Minimal tissue reaction.
- Breaks down by hydrolysis.
- Predictable absorption rate.
- Monofilament with microscopic uniform diameter.
- Soft passage through the tissues without sawing, tissue drag and trauma.
- Absence of capillarity (no support of bacterial growth).
- Excellent knot tie down.
- Pliability and softness.
- Good handling properties
- Easy and secure knotting.
- Hermetically sealed packing.
-

Guaranteed seal and product sterility.

- Indications: General soft tissue closing especially in General surgery, Gastrointestinal surgery, Ophthalmic surgery, Orthopaedics, Abdominal wall closure, Gynaecology, Plastic surgery and Urology.
- Contraindications: Cardiovascular surgery, Neurosurgery, Prosthetic devices. Being absorbable not to use when wound support for more than 6 weeks is required. Caution in patients with delayed wound healing (e.g. compromised immune system).
- Is available in a range of gauge sizes and lengths, non-needed or attached to atraumatic stainless steel needles of varying types.

| Code     | Description   | Selling unit |
|----------|---|--------------|
| 08113094 | SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 2/0 75<br>CM HR26 CYLINDER NEEDLE 26 MM   | CX. 12       |
| 08113095 | SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 3/0 75<br>CM DS24 LANCEOLATE NEEDLE 24 MM | CX. 12       |
| 08113129 | SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 0 90<br>CM HR48 CYLINDER NEEDLE 1/2 48 MM | CX. 12       |
| 08113158 | SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 2/0 75<br>CM DS24 LANCEOLATE NEEDLE 24 MM | CX. 12       |
| 08113159 | SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 0 75<br>CM DS24 LANCEOLATE NEEDLE 24 MM   | CX. 12       |
| 08113160 | SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLET 1 75<br>CM HS36 LANCEOLATE NEEDLE 36 MM   | CX. 12       |
| 08113191 | SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLETA 0<br>90CM HR36 AG/CILINDRICA 1/2 36 MM   | CX. 12       |
| 08113190 | SURGICRYL MONOFILAMENT POLYDIOXANONE VIOLETA 1/0<br>90CM HR36 AG/CILINDRICA 1/2 36 MM | CX. 12       |