

## SG - 4



### MAIN SPECIFICATIONS

- 4 compartment weekly pill box, braille marked. (LARGE).
- Helps to organize a week's medication.
- Prevents dosage duplication. Available in red and blue.
- Individual compartments hold medication for each day of the week.
- Ideal for vitamins and other medications. To open, gently push up on lid tab.
- Materials: Plastic PP.
- Sizes: Width: 216 mm Depth: 160 mm Height: 33 mm

Code	Description	Selling unit
--	SG - 4 BLUE	CX. 32
-	SG - 4 RED	CX. 32