

CUTIMED CAVITY FOAM DRESSING



MAIN SPECIFICATIONS

- The Cutimed Cavity lying between the more absorbent foam cavity wound dressings, and its increase in volume between lower, thereby protecting the wound bed. No decomposition when it gets saturated and allows for full removal, leaving no residue in the wound.
- Indications: wounds and deep cavity, for example, in venous ulcers, diabetic ulcers, pressure ulcers, wounds after excision of fistulae and abscesses, dehiscent post-surgical wounds and traumatic wounds.
- To cavity wounds with exudate levels of moderate to high.
- High absorption capacity with an increase of moderate size.
- Soft, flexible and highly malleable.
- Not adhesive.
- Microporous structure, reduces the risk of tissue growth therein and ensures a painless and non-traumatic removal.
- It can be hacked.

Code	Description	Selling unit
72621-03	CUTIMED CAVITY FOAM DRESSING 15X15 CM	CX.5
72621-01	CUTIMED CAVITY FOAM DRESSING 10X10 CM	CX. 10

FACTORMED
Produtos Médicos & Hospitalares